

XingYing\_SDK is a communication plug-in written in C++ language, in the form of a dynamic library DLL plug-in. By including the SDK into the C++ project and calling the SDK's public interface, Marker point, rigid body information and force plate including other simulation channel information can be obtained. Use reference SampleClient.

### 1. Establish connection and communication:

- (1) Create a NokovSDKClient object and enter the IP address of the XingYing target. If running on the same computer as NokovSDKClient, enter "10.1.1.198". The C++ code looks like this:

```
// The first parameter is the SDK broadcast IP of XINGYING software, see: Setting  
=>NetWork =>LocalAddress  
// ServerIP:10.1.1.198  
  
NokovSDKClient* theClient = new NokovSDKClient();  
theClient->Initialize("10.1.1.198");
```

### 2. Received data:

- (1) Set the callback function to receive the transferred data. The C++ code is as follows:  
void DataHandler(sFrameOfMocapData\* data, void\* pUserData);  
theClient->SetDataCallback(DataHandler, theClient);

### 3. Data parsing (reading Marker data and rigid-body data) :

- (1) Get all MarkerSet number:

```
int markerSetNumber = data.nMarkerSets;
```

- (2) Obtain the position information of the j-th marker of the i-th MarkerSet:

```
float point = data.MocapData[i].Markers[j][0]; // 0-x 1-y 2-z
```

- (3) Get all rigid bodies number:

```
int number = data.nRigidBodies;
```

- (4) Get the attitude information of the i-th

rigidbody(position, rotation) :

```
float point = data.RigidBodies[i].x; // x y z qx qz qw
```

- (5) Get all NamedMarker number:

```
int number = data.nLabeledMarkers;
```

- (6) Obtain the location information of the i-th NamedMarker:

```
float point = data.LabeledMarkers[i].x; // x y z
```

- (7) Get all UnNamedMarker number:

```
int number = data.nOtherMarkers;
```

- (8) Obtain the location information of the i-th UnNamedMarker:

```
float point = data.OtherMarkers[i][0]; // 0-x 1-y 2-z
```

- (9) Get the number of analog channels:

```
int number = data.nAnalogdatas;
```

(10) Get the data of the i-th analog channel:

```
float point = data.Analogdata[i];
```

(11) Get the absolute timestamp of the frame:

```
long long timestamp = data.iTimeStamp //Number of milliseconds since 1970-01-01
```

#### 4. Use of force plate:

(1) Initialize the force plate:

```
// Wait for the force plate to be initialized,0ms indicates infinite wait.
```

```
retCode = theClient->WaitForForcePlateInit(0);
```

(2) Set the callback of force plate:

```
void ForcePlateHandler(sForcePlates* pForcePlate, void* pUserData);
```

```
theClient->SetForcePlateCallback(ForcePlateHandler, theClient);
```

(3) Analyzing the data of the force plate, the force (N) on the X-axis of the i-th force table is as follows:

```
pForcePlate->ForcePlates[i].Fxyz[0];
```

#### 5. Close the connection:

```
theClient->Uninitialize();
```